Maxi Ringlight



Annular ringlight Maxi 150 mm

Features

- Continuous ring of illumination fiber within housing.
- Housing made of rugged, black anodized aluminum.
- Fiber bundle protected with flexible metal PVC sheathing.
- Large diameter ringlight for mounting to objective 58 mm outer diameter with adaptor (part no. 157 445).

Absolute shadow-free illumination for large area targets (i.e. SMD).


| Cold light source** | Part No. | Туре | Length | Inner Ø | Outer Ø | Height | Minimum working distance* | Active Ø branch | Outer Ø Jacketing | Light source adapter Ø |
|--|-------------|---------------------|------------------|-----------------|-----------------|----------------|---------------------------------|-----------------|----------------------|---------------------------------|
| | | | (L) | (ID) | (OD) | (H) | (A) | (b) | (D) | (F) |
| KL 1500 HAL KL 1600 LED KL 2500 LED KL 2500 LCD | 157 440 | Slit 150 mm Ø | 750 mm 29.53" | 150 mm 5.91" | 184 mm 7.24" | 20 mm 0.79" | 40 mm 1.57" | 9 mm 0.35" | 17 mm 0.67" | 10 mm 0.67" |


| Objective adaptor | Part No. | | |
|-------------------|----------|--|--|
| Ø 58 mm | 157 445 | | |



^{*} Ringlight usable in any working distance above given minimum values.
** and discontinued KL 1500 products (KL 1500 compact, KL 1500 LCD, KL 1500 LED, KL 1500 LED Plus)

Maxi Ringlight





Lighting and Imaging
SCHOTT AG
Hattenbergstrasse 10
55122 Mainz
Germany
Phone: +49 (0) 6131/66-7796
Fax: +49 (0) 6131/66-7850
lightingimaging@schott.com
www.schott.com/lightingimaging

